

Escape-Mattress® Economy Line

Evacuation ... in a lying position

What's an Escape-Mattress® Economy Line?

The Escape-Mattress® Economy Line is designed for people who find it difficult to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating down a spiral staircase). Evacuation mattresses are ideal for Medical Centre's and Hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients.

How does an Escape-Mattress® Economy Line work?

Place the evacuee on the Escape-Mattress® Economy with the head on to the head rest and the feet in the foot cover. After closing the safety straps the evacuee can be towed to a safe place by using the extra-long towing straps (handling min. 2 persons).

Advantages Escape-Mattress® Economy Line:

- Extra-long towing straps to help operators with improved handling practices.
- Foot cover fitted with zipper, for easier handling legs evacuee.
- Adjustable head rest.
- Compact and easy to store.
- Carrying people is no longer necessary.
- Ready for use within seconds.
- Safe evacuation down the stairs.
- Speed of descent is determined by the user.
- Easily maneuvering into the evacuation flow.
- Fits into your emergency evacuation plan.
- CE marked.

Our advice:

Make sure you have an Escape-Mattress® Economy Line available in your Organization, so you can help people who rely on you to be rescued safely. Even when you need to tackle spiral staircases.




ESCAPE MATTRESS
ECONOMY LINE

ESCAPE
MOBILITY COMPANY





Specifications	Escape-Mattress® Economy Line 
Article number:	21040
Dimensions (length x width x thick):	± 200 x 65 x 5 cm
Maximum payload:	150 kg
Weight:	± 4.8 kg incl. storage bag
Equipped with foot cover and adjustable head rest:	✓
Demonstration and advise provided:	✓
Ready for use within seconds:	✓
Speed of descent to be determined by the user:	✓
Can be folded in half:	✓
Including storage bag:	✓
CE mark:	✓
Warranty:	2 years



*Interim changes to technical specifications and printing errors reserved.
Produced by Escape Mobility Production & Innovation: Tüv approved organization.*

