

# Escape-Mattress® SOLO / FOLD / COMPACT

## Evacuation ... in a lying position

### What`s an Escape-Mattress®?

Escape-Mattress® models are designed for people who find it difficult to comfortably get into an evacuation chair, or where evacuation chairs may struggle (for example when evacuating down a spiral staircase). Evacuation mattresses are ideal for Medical Centre's and Hospitals, as they cater for complex healthcare needs. E.g. bed-bound and unconscious patients. The Escape-Mattress® models are designed for different situations to meet the needs of certain environments.

### How does an Escape-Mattress® work?

Escape-Mattress® models can be used to move evacuees along the floor (min. 2 persons) and have handles strategically placed on the side to allow easy maneuverability.

### Advantages Escape-Mattress®:

- Made of high quality materials (fire resistance classification M2).
- Extra-long towing straps to help operators with improved handling practices.
- Handgrips to assist with sideways manoeuvrability.
- Carrying people is no longer necessary.
- Ready for use within seconds.
- Safe evacuation down the stairs.
- Speed of descent is determined by the user.
- Easily maneuvering into the evacuation flow.
- Fits into your emergency evacuation plan.
- CE marked.

### Our advice:

Make sure you have an Escape-Mattress® available in your Organization, so you can help people who rely on you to be rescued safely. Even when you need to tackle spiral staircases.



Model Solo




Model Fold



Model Compact



| Specifications                                   | Escape-Mattress® Solo      | Escape-Mattress® Fold | Escape-Mattress® Compact   |
|--|----------------------------|-----------------------|----------------------------|
| Article number:                                  | 21020                      | 21010                 | 21030                      |
| Dimensions (length x width x thick):             | ± 200 x 68 x 8 cm          | ± 200 x 68 x 8 cm     | ± 200 x 68 x 6 cm          |
| Maximum payload:                                 | 150 kg                     | 150 kg                | 150 kg                     |
| Weight:  | ± 9.5 kg incl. storage bag | ± 7.2 kg              | ± 7.5 kg incl. storage bag |
| Demonstration and advise provided:               | ✓                          | ✓                     | ✓                          |
| Ready for use within seconds:                    | ✓                          | ✓                     | ✓                          |
| Speed of descent to be determined by the user:   | ✓                          | ✓                     | ✓                          |
| Fire resistance classification M2:               | ✓                          | ✓                     | ✓                          |
| Equipped with headrest, body-covers & foot sack: | ✓                          | ✓                     | ✓                          |
| Also available in narrower size:                 | ✓                          | ✓                     | -                          |
| Can be rolled up:                                | -                          | -                     | ✓                          |
| Can be folded in half:                           | -                          | ✓                     | -                          |
| Including storage bag:                           | ✓                          | Optional              | ✓                          |
| Warranty:  | 2 years                    | 2 years               | 2 years                    |



Interim changes to technical specifications and printing errors reserved.  
Produced by Escape Mobility Production & Innovation: TÜV approved organization.

